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*Place the cushion under your neck and move your head gently and evenly from side to side. With the nubbed side facing upwards, you can increase the massage effect and blood flow to the neck muscles. Place the back friend under your head with the smooth side facing upwards and stay in the resting position to relax.*

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*Place the cushion under your lumbar or thoracic spine and stretch your arms out long behind you. Alternately move one arm upwards and then turn to the side. Your gaze follows your hand.*

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*Lie on your side and lie in a long arch over your back friend. Feel the wonderful stretch in your flank. To reinforce the stretch, you can move your arm over your head to the other arm. Switch sides after approx. 30 to 60 seconds.*

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*Position both feet. Then place the back friend under your buttocks or lower back. You can tilt your pelvis slightly to the left and right from this position. Place the cushion under your thoracic spine and apply gentle pressure to the cushion.*

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### General tips for the optimal use of your Rückenfreund:


- Create a pleasant, calm and relaxed environment.
- Take your time!
- Stay in each position for at least 30 seconds to a minute.
- Breathe calmly and let your back friend carry you.

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Rückenfreund is delivered optimally inflated. If you like, you can adjust the amount of air using the enclosed needle valve pump. Please pull the valve up between your thumb and index finger as shown in the illustration. This is the easiest way to insert the needle of the pump carefully and straight into the valve. It is best to moisten the needle slightly beforehand. By the way: if you want to deflate, simply unscrew the needle from the pump, insert the needle into the valve and press lightly on the cushion. This will help you find the optimum shape for your needs.

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**EN** The instructions for use (E-IFU) can be found at [www.togu.de](http://www.togu.de) at the product.

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